

Product summary: Adult products



Motivation Active Folding Wheelchair (WM4X)



Type: Active style folding wheelchair suitable for semi-urban environments

Seat width: 360, 400, 440, 480 mm

Seat depth: 400 mm (adjustable to 360mm)

Cushion: Motivation Pressure Relief Cushion (fabric seat)

Rear wheel: 26" x 1 3/8", pneumatic tyre

Upholstery: Tension adjustable fabric and strap back and seat

Adjustments: Six standard adjustments;

- seat width
- seat depth
- backrest height (low & high)
- rear wheel balance position (active and safe)
- footrest height
- footrest angle (4° or 12°)

Net weight: 23kg

Boxed weight: 25kg

Box dimensions: 650 x 650 x 400 mm

Note: The WM4X can be re-assembled by a trained technician into any of the four seat widths and two seat depths

Motivation Rough Terrain Wheelchair (WM3)



Type: Active, three wheel wheelchair suitable for areas of rough terrain

Seat width: 360, 400, 440, 480 mm

Seat depth: 480 mm (adjustable to 400, 440 mm)

Cushion: Motivation Pressure Relief Cushion (solid seat)

Rear wheel: 26" x 1 3/8", pneumatic tyre

Backrest: Solid backrest with backrest cushion and rear pelvic pad

Adjustments:

- seat depth
- rear wheel balance position (active and safe)
- backrest angle (upright, middle and back)
- backrest height (low & high)
- footrest position (forwards & backwards)
- footrest height
- footrest angle
- rear pelvic pad height

Net weight: 22kg

Boxed weight: 25kg

Box dimensions: 890 x 645 x 275 mm

Whirlwind Rough Rider Wheelchair (RR)



Type: Active, long wheelbase, four wheel folding wheelchair suitable for rugged terrain

Seat width: 12.5" (320mm), 14" (355mm), 15.5" (395mm), 17" (430mm), 18.5" (470mm)

Seat Depth: 14" (355mm), 16" (410mm)

Each chair also comes with a 2" (5cm) extension piece

Cushion: The Comfort Cushion

Upholstery: Fabric and strap tension adjustable backrest

Rear wheel: 24" x 1.75", pneumatic tyre

Adjustments:

- rear wheel balance position (5)
- backrest height (low, middle and high)
- seat depth (utilising a 50mm seat depth extension)
- footrest height

Net weight: 21.5kg

Boxed weight: 25kg (approximately)

Boxed dimensions: 12 x 36 x 26" (310 x 920 x 665 mm)

Note: A Moderate Pressure Management Cushion Lift can be ordered with the Comfort Cushion



LDSC Standard wheelchair



Type: Standard wheelchair, folding with removable footplates

Seat width: 300mm, 350mm, 410mm, 450mm, 500mm, 550mm, 600mm

Seat depth: 330 – 430mm adjustable

Cushion: Basic foam cushion

Rear wheel: 24" x 1 3/8", solid polyurethane tyre

Upholstery: Nylon, standard upholstery

Adjustments:

- backrest height
- footrest height
- arm rest height
- floor to seat height

Box dimensions: 710 x 400 x 660 mm

LDSC Active wheelchair



Type: Active wheelchair, seat back folds down

Seat width: 300 mm, 350 mm, 400 mm,
450 mm, 500 mm, 550 mm

Seat depth: 330 – 430mm adjustable

Cushion: Basic foam cushion

Rear wheel: 24" x 1 3/8", solid pneumatic tyre, quick release axles

Upholstery: Nylon

Adjustments:

- rear wheel position
- backrest height, angle
- footrest height, angle
- arm rest height, and can be fixed, removable or swing away
- floor to seat height
- optional anti-tipper

Box dimensions: 650 x 750 x 490 mm

LDSC All-terrain wheelchair



Type: Rough terrain wheelchair, seat back folds down

Seat width: 300 mm, 350 mm, 400 mm,
450 mm, 500 mm, 550 mm

Seat depth: 330 – 430mm adjustable

Cushion: Basic foam cushion

Rear wheel: 24" x 1 3/8", solid pneumatic tyre, quick release axles

Upholstery: Nylon

Adjustments:

- rear wheel position
- backrest height, angle
- footrest height, angle
- arm rest height, and can be fixed, removable or swing away
- floor to seat height
- optional anti-tipper

Motivation Clip on Trike (TRCO)



Type: The Clip-on Tricycle Attachment is a quick release clip-on, clip-off attachment, which can be fitted to the Motivation Rough Terrain (WM3) wheelchair. It converts the wheelchair into a tricycle allowing the user to travel greater distances more effectively

Front wheel: 20" x 1.75", pneumatic tyre

Net weight: 13kg (plus 22kg for the WM3)

Boxed weight: 16kg

Box dimensions: 1325 x 195 x 550 mm



Product summary: Cushions

Motivation Pressure Relief Cushion (fabric seat)



A moulded contoured foam cushion with medium pressure relief and positioning.

Each cushion can be cut to size and has a size-adjustable, water resistant, nylon cover.

Uncut size: 480 x 440 mm

Net weight: 1.3kg (uncut)

Note: Cushions sold in boxes of 10 x units, with or without cushion covers.



Motivation Pressure Relief Cushion (solid seat)



A moulded contoured foam cushion with medium pressure relief and positioning.

Each cushion can be cut to size and has a size-adjustable, water resistant, nylon cover.

Uncut size: 480 x 440 mm

Net weight: 1.3kg (uncut)

Note: Cushions sold in boxes of 10 x units, with or without cushion covers



Jarik Fluid Wheelchair Cushion



A combination foam and fluid cushion, providing moderate to high level of pressure relief. The cushion has a contoured foam base with a fluid pad sitting in the rear pelvic well. There is a waterproof base cover, with a stretchy top cover, which can be easily removed for washing.

Cushion sizes: 14" x 17", 14" x 19", 16" x 17", 16" x 19", 18" x 17", 18" x 19", 20" x 19".